**REGULATORY:** Class I device. Complies with MDR (EU) 2017/745 and manufactured according to the procedures of BIONEN Quality System certified ISO 13485.

**GUARANTEE AND LIMITATION OF LIABILITY:**

BIONEN guarantees that the device is free of material and manufacturing defects for a period of twelve months from the date of purchase. The sole responsibility of BIONEN and our dealers is limited to the replacement or repair of the product, in the opinion of BIONEN, without charge for parts or work, if it is proven that one of the parts is defective in manufacturing, performance or materials, during warranty period. Under no circumstances will BIONEN or its dealers be liable for any loss of earnings or for damages, direct, indirect or accidental, including loss of profit, damage to property, or personal injury arising from the use, or from the inability to use the product. This warranty only covers the original purchaser and supersedes all other warranties or prior agreements, either expressed or implied. This warranty is void if the product is used for purposes other than those for which it is intended, or if it is subject to abuse, misuse, tampering, negligence, or unauthorized changes. Use of this product constitutes acceptance of the full warranty. Notify the BIONEN Quality Assurance Service for any malfunction or defect found on this product.

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## MUSCLE STIMULATOR

| REF.                 | VERSIONS                                  |
|----------------------|-------------------------------------------|
| 0038.0xx             | Flat bridge stimulator with belt / handle |
| 0038.02x<br>0038.03x | Flat bridge stimulator with belt and pads |
| 0038.2xx             | Stimulator with velcro fastening          |
| 0039.0xx             | Bipolar muscle stimulator                 |
| 0039.20x             | Monopolar muscle stimulator               |

**INTENDED USE:** EMG - EP

**INSTRUCTION FOR USE:** Lean the stimulator's electrode on the point to be stimulated with a light pressure.

In case of stimulators with pads, dip the pads in a 0,9% physiological saline solution. Clean up patient's skin and apply the stimulator using the belt.

Connect the plugs to the device, making sure to reach the best contact impedance.

**CLEANING INSTRUCTIONS:** Carefully clean the stimulator after each use with a wet cloth/gauze. Do not dip into any liquid.

In case of stimulators with pads, after each test, clean the felts with water and neutral detergent.

**DISINFECTION:** Disinfect using either Amocid, Lysoformin, Alhyndex, IndicinGG or equivalent products, at high concentration too, without immersing the device into the liquids.

**STERILIZATION METHOD:** NON-STERILE device that does not need to be sterilized before use.

**LIFETIME:** After 5 years from sale, is not guaranteed the proper functionality of the device.

**STORAGE:** Particular precautions for storing are not necessary but it's best to put the electrodes in its box after washing

**CONTRAINDICATIONS AND INTERACTIONS:** none

**POSSIBLE TOXICITY DECLARED:** none. LATEX FREE

**WARNINGS AND PRECAUTIONS:** For professional use only.

**TRANSPORT AND DISPOSAL**

There are no special modes of transport and disposal.

**HANDLING AND STORAGE:** Do not overlap excessive loads.

**WASTE DISPOSAL:** manage the product as sanitary waste in conformity with the regulations in force